

Silonite™ Canisters with the **NEW** TrueSeal™ Valve



15L Canister
w/ TSG-01

6L Canister
w/ TSG-01

3.2L Canister
w/ TSG-01

2.7L Canister
w/ TSG-01

2L Canister
w/ TSG-01

1.4L Canister
w/ TSG-01

1L Canister
w/ TSG-01

0.01–100 PPB Trace Level

Entech's Silonite™ canisters feature a large volume capacity for detection of volatile chemicals down to the low part per trillion range. An inert and durable internal Silonite™ coating provides a high-quality, long-term sample storage solution. Losses in the valve are avoided by using our new TrueSeal™ Valve which includes Entech's new "sure-seal" technology. An integrated valve guard is securely welded to the canister for superior light-weight valve protection without any stresses to the valve stem associated with heavy "strapped-type" valve guards. These canisters are certified to meet or exceed the technical specifications required for EPA methods TO-14a and TO-15. High quality performance is verified for EVERY canister with our demanding chemical inertness tests.

Features

- ◆ **Silonite™ Surface Coating**
We believe that each and every canister should be as inert and corrosion resistant as possible. That is why all of our canisters feature Silonite™.
- ◆ **Welded Valve Guard**
Unlike "strap" designs, our TrueSeal™ guard provides superior valve protection and far less added weight.
- ◆ **Extended Range Canister Sampling**
Our Silonite™ Canisters are more inert than SUMMA® or electropolished canisters, allowing recovery of a greater range of chemicals – including reduced sulfur compounds.

The 3.2L and 2.7L are lighter alternatives for canister sampling featuring minimal weight and maximum portability. 6L and 15L Canisters are ideal for large volume sampling and lab standards.

Description	Unit	Part #
1L Can Silonite™ TrueSeal™ Valve	EA	29-MC10LSV
1L Can Silonite™ TrueSeal™ & Compound® Gauge	EA	29-MC10LSVG
1.4L Can Silonite™ TrueSeal™ Valve	EA	29-MC14LSV
1.4L Can Silonite™ TrueSeal™ & Compound® Gauge	EA	29-MC14LSVG
2L Can Silonite™ TrueSeal™ Valve	EA	29-MC2LSV
2L Can Silonite™ TrueSeal™ & Compound® Gauge	EA	29-MC2LSVG
2.7L Can Silonite™ TrueSeal™ Valve	EA	29-10272A
2.7L Can Silonite™ TrueSeal™ & Compound® Gauge	EA	29-10272VG
3.2L Can Silonite™ TrueSeal™ Valve	EA	29-10322
3.2L Can Silonite™ TrueSeal™ & Compound® Gauge	EA	29-10322VG
6L Can Silonite™ TrueSeal™ Valve	EA	29-10622
6L Can Silonite™ TrueSeal™ & Compound® Gauge	EA	29-10622VG
15L Can Silonite™ TrueSeal™ Valve	EA	29-11522
15L Can Silonite™ TrueSeal™ & Compound® Gauge	EA	29-11522VG