

## Silonite-XL™ Canisters with the **NEW** TrueSeal™ Valve



### 0.01–100 PPB Trace Level

Entech's Silonite-XL™ canisters feature a large volume capacity for detection of volatile chemicals down to the low part per trillion range. An inert and durable internal Silonite-XL™ coating provides a high-quality, long-term sample storage solution. Losses in the valve are avoided by using our new TrueSeal™ Valve which includes Entech's new "sure-seal" technology. An integrated valve guard is securely welded to the canister for superior light-weight valve protection without any stresses to the valve stem associated with heavy "strapped-type" valve guards. These canisters are certified to meet or exceed the technical specifications required for EPA methods TO-14a and TO-15. High quality performance is verified for EVERY canister with our demanding chemical inertness tests.

### Features

- Silonite-XL™ Surface Coating**  
*We believe that each and every canister should be as inert and corrosion resistant as possible. That is why all of our canisters feature Silonite-XL™.*
- Welded Valve Guard**  
*Unlike "strap" designs, our TrueSeal™ guard provides superior valve protection and far less added weight.*
- Extended Range Canister Sampling**  
*Our Silonite-XL™ Canisters are more inert than SUMMA® or electropolished canisters, allowing recovery of a greater range of chemicals – including reduced sulfur compounds.*

*The 3.2L and 2.7L are lighter alternatives for canister sampling featuring minimal weight and maximum portability. 6L and 15L Canisters are ideal for large volume sampling and lab standards.*

Description	Unit	Part #
1L Can Silonite-XL™ TrueSeal™ Valve	EA	29-MC10LSV
1L Can Silonite-XL™ TrueSeal™ & Compound® Gauge	EA	29-MC10LSVG
1.4L Can Silonite-XL™ TrueSeal™ Valve	EA	29-MC14LSV
1.4L Can Silonite-XL™ TrueSeal™ & Compound® Gauge	EA	29-MC14LSVG
2.5L Can Silonite-XL™ TrueSeal™ Valve	EA	29-10252
2.5L Can Silonite-XL™ TrueSeal™ & Compound® Gauge	EA	29-10252VG
2.7L Can Silonite-XL™ TrueSeal™ Valve	EA	29-10272A
2.7L Can Silonite-XL™ TrueSeal™ & Compound® Gauge	EA	29-10272VG
3.2L Can Silonite-XL™ TrueSeal™ Valve	EA	29-10322
3.2L Can Silonite-XL™ TrueSeal™ & Compound® Gauge	EA	29-10322VG
6L Can Silonite-XL™ TrueSeal™ Valve	EA	29-10622
6L Can Silonite-XL™ TrueSeal™ & Compound® Gauge	EA	29-10622VG
15L Can Silonite-XL™ TrueSeal™ Valve	EA	29-11522
15L Can Silonite-XL™ TrueSeal™ & Compound® Gauge	EA	29-11522VG